

Anti-Crybb1 Rabbit Polyclonal Antibody For Wb, Elisa - P53674

Item Number: CM0025149



Introduction

Rabbit polyclonal antibody targeting human CRYBB1 (beta-crystallin B1), validated for WB and ELISA. Suitable for eye lens protein studies and cataract disease variant analysis. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	CRYBB1 Rabbit pAb
Remarks/Alias	CATCN3; CTRCT17; CRYBB1
Species	Human
Gene ID (Human)	1414
Gene ID	1414
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human CRYBB1 (NP_001878.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P53674
Protein Weight	28kDa
Shipping	Ice bag