


Anti-Crybb2 Rabbit Polyclonal Antibody For Wb, Elisa - P43320

Item Number: CM0026541



Introduction

Rabbit polyclonal antibody targeting beta-crystallin B2 (CRYBB2). Validated for Western blot and ELISA applications, detecting endogenous CRYBB2 in human, mouse, and rat samples. Molecular weight ~23 kDa.

[Learn More](#)

Parameter	Value
Product Name	CRYBB2 Rabbit pAb
Remarks/Alias	CCA2, CRYB2, CRYB2A, CTRCT3, D22S665, CRYBB2
Species	Human
GeneID (Human)	1415
GeneID	1415
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human CRYBB2 (NP_000487.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43320
Protein Weight	23kDa
Shipping	Ice bag