

Anti-Crygc Polyclonal Antibody For Wb, If/Icc, Elisa - P07315

Item Number: CM0030317



Introduction

Rabbit polyclonal antibody targeting human gamma-crystallin C (CRYGC). Validated for Western blot, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Useful for eye lens protein and cataract research.

[Learn More](#)

Parameter	Value
Product Name	CRYGC Rabbit pAb
Remarks/Alias	CCL; CRYG3; CTRCT2; CRYGC
Species	Human
GeneID (human)	1420
GeneID	1420
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07315
Protein Weight	21kDa
Shipping	Ice bag