


Anti-Crygs Polyclonal Antibody For Wb, Elisa - P22914

Item Number: CM0026749



Introduction

High-quality rabbit polyclonal antibody targeting human CRYGS (gamma-crystallin S). Suitable for WB and ELISA applications, with demonstrated cross-reactivity in human, mouse, and rat samples. Ideal for lens protein and cataract research.

[Learn More](#)

Parameter	Value
Product Name	CRYGS Rabbit pAb
Remarks/Alias	CRYG8; CTRCT20; CRYGS
Species	Human
GeneID (Human)	1427
GeneID	1427
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human CRYGS (NP_060011.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22914
Protein Weight	21kDa
Shipping	Ice bag