

Anti-Crym Polyclonal Antibody For Wb, Elisa - Q14894

Item Number: CM0030903



Introduction

Rabbit polyclonal antibody targeting CRYM (ketimine reductase mu-crystallin), validated for WB and ELISA. Reacts with human, mouse, and rat samples. Suitable for metabolism and deafness research.

[Learn More](#)

Parameter	Value
Product Name	CRYM Rabbit pAb
Remarks/Alias	THBP; DFNA40; CRYM
Species	Human
GeneID (Human)	1428
GeneID	1428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-314 of human CRYM (NP_001879.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14894
Protein Weight	34kDa
Shipping	Ice bag