

Anti-Cryz Polyclonal Antibody For Wb, Elisa - Q08257

Item Number: CM0019456



Introduction

Rabbit polyclonal antibody against CRYZ (zeta-crystallin), KO validated. Suitable for Western blot and ELISA; cross-reacts with human, mouse, and rat. Useful for mRNA stabilization, quinone reductase, and detoxification research. SWISS Q08257.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CRYZ Rabbit pAb
Remarks/Alias	CRYZ; YZ
Species	Human
Gene ID (Human)	1429
Gene ID	1429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-329 of human CRYZ (NP_001123514.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08257
Protein Weight	35kDa
Shipping	Ice bag