


Anti-Csf1 Rabbit Polyclonal Antibody For Wb And Elisa - P09603

Item Number: CM0020321



Introduction

This rabbit polyclonal antibody against human CSF1 (M-CSF) is validated for WB and ELISA. It cross-reacts with mouse and rat, and detects a ~60 kDa protein. Useful for studies of macrophage differentiation and immune response.

[Learn More](#)

Parameter	Value
Product Name	CSF1 Rabbit pAb
Remarks/Alias	MCSF; CSF-1; CSF1
Species	Human
Gene ID (Human)	1435
Gene ID	1435
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-190 of human CSF1 (NP_757351.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09603
Protein Weight	60kDa
Shipping	Ice bag