


Anti-Csf1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P09603

Item Number: CM0015059



Introduction

Anti-CSF1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA with human samples. Essential for macrophage colony-stimulating factor 1 studies in hematopoietic differentiation, innate immunity, and bone development.

[Learn More](#)

Parameter	Value
Product Name	CSF1 Rabbit pAb
Remarks/Alias	MCSF; CSF-1; CSF1
Species	Human
GeneID	1435
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-190 of human CSF1 (NP_757351.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P09603
Protein Weight	60kDa
Shipping	Ice bag