

Anti-Csf2Ra Rabbit Polyclonal Antibody For If/Icc, Elisa - P15509

Item Number: CM0027269



Introduction

Rabbit polyclonal antibody targeting human CSF2RA (GM-CSF receptor alpha). Validated for IF/ICC and ELISA; cross-reacts with human and mouse. Ideal for hematopoietic cell and receptor studies.

[Learn More](#)

Parameter	Value
Product Name	CSF2RA Rabbit pAb
Remarks/Alias	GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; alphaGMR; GMR-alpha; GMCSFR-alpha; GM-CSF-R-alpha; CSF2RA
Species	Human
GeneID (Human)	1438
GeneID	1438
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-150 of human CSF2RA (NP_001155001.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P15509
Protein Weight	46kDa
Shipping	Ice bag