

Anti-Csf2Rb Rabbit Polyclonal Antibody For Wb, Elisa - P32927

Item Number: CM0013860



Introduction

Rabbit polyclonal antibody against human CSF2RB (CD131), suitable for western blot (WB) and ELISA. Cross-reacts with human and mouse. Facilitates research on IL-3, IL-5, and GM-CSF signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	CSF2RB Rabbit pAb
Remarks/Alias	CD131; IL3RB; IL5RB; SMDP5; CDw131; betaGMR; CSF2RB
Species	Human
GeneID (Human)	1439
GeneID	1439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-443 of human CSF2RB (NP_000386.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32927
Protein Weight	97kDa
Shipping	Ice bag