

Anti-Csf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P09919

Item Number: CM0024984



Introduction

CSF3 Rabbit Polyclonal Antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat CSF3; 22 kDa. Recombinant immunogen; suitable for granulocyte colony-stimulating factor studies in hematopoiesis and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CSF3 Rabbit pAb
Remarks/Alias	GCSF; CSF3OS; C17orf33; CSF3
Species	Human
GeneID (human)	1440
GeneID	1440
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of human CSF3 (NP_000750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09919
Protein Weight	22kDa
Shipping	Ice bag