


Anti-Csf3R Polyclonal Antibody For Wb, Elisa - Q99062

Item Number: CM0028870



Introduction

Anti-CSF3R rabbit polyclonal antibody for Western blot and ELISA, specifically detecting human granulocyte colony-stimulating factor receptor. Ideal for studying neutrophilic lineage, granulocytic maturation, and leukemia models.

[Learn More](#)

Parameter	Value
Product Name	CSF3R Rabbit pAb
Remarks/Alias	SCN7; CD114; GCSFR; CSF3R
Species	Human
GeneID (human)	1441
GeneID	1441
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 650-750 of human CSF3R (NP_000751.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99062
Protein Weight	92kDa
Shipping	Ice bag