


Anti-Csk (C-Src Kinase) Rabbit Polyclonal Antibody For Wb, Elisa - P41240

Item Number: CM0018211



Introduction

CSK Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, rat CSK, a key Src-family kinase regulator. Immunogen: recombinant human CSK (aa 1-200). Ideal for signaling and immunity research.

[Learn More](#)

Parameter	Specification
Product Name	CSK Rabbit pAb
Remarks/Alias	CSK
Species	Human
GeneID (Human)	1445
GeneID	1445
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CSK (NP_001120662.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41240
Protein Weight	51kDa
Shipping	Ice bag