


Anti-Csk Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P41240

Item Number: CM0002557



Introduction

High-quality rabbit monoclonal antibody against human CSK, validated for Western blot, immunohistochemistry (paraffin), and ELISA. Detects CSK in human, mouse, and rat samples, supporting research into Src kinase regulation and immune signaling.

[Learn More](#)

Parameter	Value
Product Name	CSK Rabbit mAb
Remarks/Alias	CSK; tyrosine-protein kinase CSK
Species	Human
Gene ID (Human)	1445
Gene ID	1445
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CSK (P41240).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41240
Protein Weight	51kDa
Shipping	Ice bag