


Anti-Csk Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P41240

Item Number: CM0023657



Introduction

Rabbit polyclonal antibody targeting human CSK tyrosine kinase, validated for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for investigating CSK-mediated regulation of Src-family kinases in immune response and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CSK Rabbit pAb
Remarks/Alias	CSK
Species	Human
Gene ID	1445 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 188-450 of human CSK (NP_004374.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41240
Protein Weight	51kDa
Shipping	Ice bag