

Anti-Csmd3 Polyclonal Antibody For Wb, Elisa - Q7Z407

Item Number: CM0015468

Introduction

Rabbit polyclonal antibody targeting human CSMD3, a large transmembrane protein involved in dendrite development. Validated for WB and ELISA, with reactivity confirmed in human samples. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CSMD3 Rabbit pAb
Remarks/Alias	CSMD3
Species	Human
GeneID (Human)	114788
GeneID	114788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-440 of human CSMD3 (NP_937756.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z407
Protein Weight	406kDa
Shipping	Ice bag