


Anti-Csnk1E Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49674

Item Number: CM0018770



Introduction

High-quality rabbit polyclonal antibody against human CSNK1E (Casein kinase I isoform epsilon), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for circadian clock and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CSNK1E Rabbit pAb
Remarks/Alias	CK1e; HCK1E; CK1epsilon; 1E
Species	Human
GeneID (human)	1454
GeneID	1454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 247-416 of human CSNK1E (NP_689407.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49674
Protein Weight	47kDa
Shipping	Ice bag