


Anti-Csnk1G1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hcp0

Item Number: CM0028073



Introduction

CSNK1G1 Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human casein kinase I isoform gamma-1, a serine/threonine kinase involved in Wnt signaling. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	CSNK1G1 Rabbit pAb
Remarks/Alias	CK1gamma1; CSNK1G1
Species	Human
GeneID (human)	53944
GeneID	53944
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 298-422 of human CSNK1G1 (NP_071331.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HCP0
Protein Weight	49kDa
Shipping	Ice bag