

Anti-Csnk1G3 Polyclonal Antibody For Wb, Elisa - Q9Y6M4

Item Number: CM0016586



Introduction

CamelBio's rabbit anti-human CSNK1G3 polyclonal antibody (pAb) detects casein kinase I gamma-3 via WB and ELISA, with cross-reactivity to mouse and rat, supporting Wnt signaling and kinase research.

[Learn More](#)

Parameter	Value
Product Name	CSNK1G3 Rabbit pAb
Remarks/Alias	CSNK1G3L; CKI-gamma 3; CSNK1G3
Species	Human
GeneID (human)	1456
GeneID	1456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 304-423 of human CSNK1G3 (NP_001026982.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y6M4
Protein Weight	51kDa
Shipping	Ice bag