

Anti-Csnk2A2 Polyclonal Antibody For Wb, If/Icc, Elisa - P19784

Item Number: CM0023027



Introduction

Rabbit polyclonal antibody against human CSNK2A2 (CK2 α'), validated for western blot, immunofluorescence/ICC, and ELISA. Detects human, mouse, and rat CSNK2A2. Ideal for studying cell cycle, apoptosis, and transcriptional regulation pathways.

[Learn More](#)

Parameter	Value
Product Name	CSNK2A2 Rabbit pAb
Remarks/Alias	CK2A2; CSNK2A1; CK2 α' ; CSNK2A2
Species	Human
GeneID (Human)	1459
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CSNK2A2 (NP_001887.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19784
Protein Weight	41kDa
Shipping	Ice bag