

Anti-Csnk2A2 Rabbit Polyclonal Antibody For Wb And Elisa - P19784

Item Number: CM0023459



Introduction

Rabbit polyclonal antibody targeting human and mouse CSNK2A2, suitable for Western blot and ELISA. Recombinant protein immunogen ensures specificity. Useful for research in kinase signaling, cell cycle, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	CSNK2A2 Rabbit pAb
Remarks/Alias	CK2A2; CSNK2A1; CK2alpha'; CSNK2A2
Species	Human
GeneID (human)	1459
GeneID	1459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human CSNK2A2 (NP_001887.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19784
Protein Weight	41kDa
Shipping	Ice bag