

Anti-Cspp1 Polyclonal Antibody For Wb, Elisa - Q1Msj5

Item Number: CM0028688

Introduction

Rabbit polyclonal antibody targeting human CSPP1, validated for WB and ELISA. Detects centrosome and spindle pole-associated protein 1 in human and mouse samples. Ideal for cell cycle and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	CSPP1 Rabbit pAb
Remarks/Alias	CSPP; CSPP-L; JBTS21; CSPP1
Species	Human
GeneID (Human)	79848
GeneID	79848
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 962-1221 of human CSPP1 (NP_079066.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q1MSJ5
Protein Weight	146kDa
Shipping	Ice bag