

Anti-Csrnp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96S65

Item Number: CM0024804

Introduction

Rabbit polyclonal antibody targeting human CSRNP1 (AXUD1, TAIP-3), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for apoptosis and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CSRNP1 Rabbit pAb
Remarks/Alias	AXUD1; URAX1; TAIP-3; CSRNP-1; FAM130B; CSRNP1
Species	Human
GeneID	64651
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CSRNP1 (NP_149016.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96S65
Protein Weight	64kDa
Shipping	Ice bag