

Anti-Csrp1 Polyclonal Antibody For Wb, Ip, Elisa, Rip - P21291

Item Number: CM0017940

Introduction

Rabbit polyclonal antibody against human CSRP1, suitable for WB, IP, ELISA, and RIP. Cross-reacts with mouse. Ideal for investigating CSRP1's role in neuronal development and nuclear functions.

[Learn More](#)

Parameter	Value
Product Name	CSRP1 Rabbit pAb
Remarks/Alias	CRP; CRP1; CSRP; CYRP; D1S181E; HEL-141; HEL-S-286; CSRP1
Species	Human
Gene ID (Human)	1465
Gene ID	1465
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA, RIP
Cross Reactivity	Human, Mouse
SWISS	P21291
Protein Weight	21kDa
Shipping	Ice bag