

Anti-Csrp2 Polyclonal Antibody For Wb, Elisa - Q16527

Item Number: CM0024937

Introduction

Rabbit polyclonal antibody targeting human CSRP2 (CRP2/LMO5/SmLIM), reactive with human and mouse. For WB and ELISA applications in vascular biology and smooth muscle cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	CSRP2 Rabbit pAb
Aliases	CRP2; LMO5; SmLIM; CSRP2
Species	Human
GeneID (Human)	1466
GeneID	1466
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human CSRP2 (NP_001312.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16527
Protein Weight	21kDa
Shipping	Ice bag