

Anti-Csrp2Bp Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H8E8

Item Number: CM0028707

Introduction

CSRP2BP Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA with reactivity to human and rat. Targets the Cysteine-rich protein 2-binding protein (KAT14/ATAC2), a subunit of the ATAC histone acetyltransferase complex involved in chromatin remodeling and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	CSRP2BP Rabbit pAb
Remarks/Alias	ATAC2; CRP2BP; CSRP2BP; PRO1194; dj717M23.1
Species	Human
GeneID (Human)	57325
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 523-782 of human CSRP2BP (NP_065397.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H8E8
Protein Weight	89kDa
Shipping	Ice bag