


Anti-Csrp3 Rabbit Polyclonal Antibody - P50461

Item Number: CM0028185



Introduction

Anti-CSR3 rabbit polyclonal antibody suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential for muscle differentiation and cardiac function. Ideal for cardiomyopathy research.

[Learn More](#)

Parameter	Value
Product Name	CSR3 Rabbit pAb
Remarks/Alias	CLP; MLP; CRP3; LMO4; CMD1M; CMH12; CSR3
Species	Human
GeneID (Human)	8048
GeneID	8048
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human CSR3 (NP_003467.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50461
Protein Weight	21kDa
Shipping	Ice bag