

## Anti-Cst1 Polyclonal Antibody For Wb, Elisa - P01037

Item Number: CM0030771

### Introduction

CST1 Rabbit pAb detects human Cystatin-SN (~14 kDa), a secreted cysteine protease inhibitor. Validated for WB and ELISA; suitable for studies of saliva, tears, urine, and seminal fluid.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CST1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | CST1                                                                                                             |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 1469                                                                                                             |
| Gene ID          | 1469                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-141 of human CST1 (NP_001889.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | P01037                                                                                                           |
| Protein Weight   | 14kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |