

Anti-Cst6 (Cystatin-M) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15828

Item Number: CM0027740



Introduction

High-affinity rabbit polyclonal antibody targeting human CST6 (Cystatin-M), validated for WB, IF/ICC, and ELISA. Ideal for studying epidermal cornification and hair follicle biology. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CST6 Rabbit pAb
Remarks/Alias	ECTD15; CST6
Species	Human
GeneID (human)	1474
GeneID	1474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-149 of human CST6 (NP_001314.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15828
Protein Weight	17kDa
Shipping	Ice bag