

Anti-Cst8 Rabbit Polyclonal Antibody For Wb, Elisa - O60676

Item Number: CM0017140

Introduction

Anti-CST8 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human CST8 and cross-reacts with mouse. Suitable for studying sperm development and epididymal function. SWISS: O60676.

[Learn More](#)

Parameter	Value
Product Name	CST8 Rabbit pAb
Remarks/Alias	CRES; CTES5; CST8
Species	Human
GeneID (Human)	10047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 53-142 of human CST8 (NP_005483.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O60676
Protein Weight	16kDa
Shipping	Ice bag