

Anti-Cstb Polyclonal Antibody For Wb, Elisa - P04080

Item Number: CM0027847



Introduction

Rabbit polyclonal antibody targeting human CSTB (Cystatin-B), a thiol protease inhibitor. Validated for WB and ELISA, cross-reacts with human and rat. Ideal for neurodegeneration, epilepsy, and protease inhibition research.

[Learn More](#)

Parameter	Value
Product Name	CSTB Rabbit pAb
Remarks/Alias	PME; ULD; CST6; EPM1; STFB; CPI-B; EPM1A; CSTB
Species	Human
GeneID (Human)	1476
GeneID	1476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human CSTB (NP_000091.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04080
Protein Weight	11kDa
Shipping	Ice bag