

Anti-Cstb Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04080

Item Number: CM0022576



Introduction

Rabbit polyclonal antibody targeting human cystatin-B (CSTB), an 11kDa intracellular thiol proteinase inhibitor. Suitable for WB, IHC-P, ELISA. Immunogen: recombinant protein (aa 1-98).

[Learn More](#)

Parameter	Value
Product Name	CSTB Rabbit pAb
Remarks/Alias	PME; ULD; CST6; EPM1; STFB; CPI-B; EPM1A; CSTB
Species	Human
Gene ID	1476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human CSTB (NP_000091.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P04080
Protein Weight	11kDa
Shipping	Ice bag