

Anti-Cstf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P33240

Item Number: CM0027721



Introduction

Rabbit polyclonal antibody targeting CSTF2 (CstF-64), a key factor in mRNA 3' end processing. Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for polyadenylation and cleavage research.

[Learn More](#)

Parameter	Value
Product Name	CSTF2 Rabbit pAb
Remarks/Alias	CstF-64; CSTF2
Species	Human
GeneID (Human)	1478
GeneID	1478
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human CSTF2 (NP_001316.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33240
Protein Weight	61kDa
Shipping	Ice bag