

## Anti-Cstf2T Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H0L4

Item Number: CM0026348

### Introduction

CSTF2T rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human and mouse CSTF2T, a nuclear protein involved in AAUAAA-independent mRNA polyadenylation in germ cells.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | CSTF2T Rabbit pAb                                                                                                   |
| Remarks/Alias    | CstF-64T; CSTF2T                                                                                                    |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 23283                                                                                                               |
| GeneID           | 23283                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 420-560 of human CSTF2T (NP_056050.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                        |
| SWISS            | Q9H0L4                                                                                                              |
| Protein Weight   | 64kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |