

Anti-Ctbp1 Monoclonal Antibody For Wb, Elisa - Q13363

Item Number: CM0003910



Introduction

Rabbit monoclonal antibody against CtBP1, validated for Western blot and ELISA. Reacts with human, mouse, rat samples. Suitable for studying transcription regulation, Golgi structure, and BAT differentiation.

[Learn More](#)

Parameter	Value
Product Name	CtBP1 Rabbit mAb
Remarks/Alias	BARS, HADDTS, CtBP1
Species	Human
GeneID (Human)	1487
GeneID	1487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 301-440 of human CtBP1 (Q13363).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13363
Protein Weight	48 kDa
Shipping	Ice bag