

Anti-Ctbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13363

Item Number: CM0020462



Introduction

High-quality rabbit polyclonal antibody targeting human CtBP1, validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for studying transcriptional corepression, Golgi structure, and BAT differentiation.

[Learn More](#)

Parameter	Specification
Product Name	CtBP1 Rabbit pAb
Remarks/Alias	BARS; HADDTS; CtBP1
Species	Human
GeneID	1487 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-250 of human CtBP1 (NP_001319.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13363
Protein Weight	48kDa
Shipping	Ice bag