


Anti-Ctbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13363

Item Number: CM0016416



Introduction

Human CtBP1 rabbit polyclonal antibody validated for WB and ELISA, reacting with human, mouse, and rat. Ideal for studying transcriptional repression, Golgi structure, and brown adipose differentiation. - Q13363

[Learn More](#)

Parameter	Value
Product Name	CtBP1 Rabbit pAb
Remarks/Alias	BARS; HADDT5; CtBP1
Species	Human
Gene ID (Human)	1487
Gene ID	1487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-353 of CtBP1 (NP_001319.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13363
Protein Weight	48kDa
Shipping	Ice bag