

## Anti-Ctbp2 Polyclonal Antibody For Wb, If/Icc, Elisa - P56545

Item Number: CM0024561



### Introduction

High-quality CTBP2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat CTBP2 (49 kDa). Ideal for studying transcriptional repression, brown adipocyte differentiation, and synaptic organization.

[Learn More](#)

| Parameter        | Specification                                                                                                              |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CTBP2 Rabbit pAb                                                                                                           |
| Remarks/Alias    | CTBP2                                                                                                                      |
| Species          | Human                                                                                                                      |
| GeneID (human)   | 1488                                                                                                                       |
| GeneID           | 1488                                                                                                                       |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 345-445 of human CTBP2 (NP_001320.1). |
| Source           | Rabbit                                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                                      |
| Application      | WB, IF/ICC, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                          |
| SWISS            | P56545                                                                                                                     |
| Protein Weight   | 49kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |