

Anti-Ctbp2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P56545

Item Number: CM0002233



Introduction

Highly specific CTBP2 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat. Ideal for transcriptional corepressor and synapse scaffold studies.

[Learn More](#)

Parameter	Value
Product Name	CTBP2 Rabbit mAb
Remarks/Alias	CTBP2; C-terminal-binding protein 2
Species	Human
GeneID (Human)	1488
GeneID	1488
Immunogen	Synthetic Peptide, A synthetic peptide corresponding to a sequence within amino acids 346-445 of human CTBP2 (P56545).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56545
Protein Weight	49kDa
Shipping	Ice bag