

Anti-Ctcf Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Chip, Elisa - P49711

Item Number: CM0015927



Introduction

High-quality rabbit anti-CTCF polyclonal antibody validated for WB, IP, IF/ICC, IHC-P, ChIP, ELISA. Reacts with human, mouse, rat. CTCF is a key chromatin insulator and transcriptional regulator. Ideal for chromatin research.

[Learn More](#)

Parameter	Specification
Product Name	CTCF Rabbit pAb
Remarks/Alias	MRD21; FAP108; CFAP108; CTCF
Species	Human
Gene ID (Human)	10664
Gene ID	10664
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CTCF (NP_006556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49711
Protein Weight	83kDa
Shipping	Ice bag