

Anti-Ctcf1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ni51

Item Number: CM0024703

Introduction

High-quality rabbit polyclonal antibody targeting human CTCFL (BORIS/CT27), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying transcriptional regulation, imprinting, and cancer-testis antigens.

[Learn More](#)

Parameter	Value
Product Name	CTCFL Rabbit pAb
Remarks/Alias	CT27; BORIS; CTCF-T; HMGB1L1; dj579F20.2; CTCFL
Species	Human
Gene ID (Human)	140690
Gene ID	140690
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CTCFL (NP_001255970.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NI51
Protein Weight	76kDa
Shipping	Ice bag