

Anti-Ctgf Polyclonal Antibody For Wb, If/Icc, Elisa - P29279

Item Number: CM0018192

Introduction

Rabbit polyclonal antibody against human CTGF (CCN2), validated for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with rat. Ideal for connective tissue, osteogenesis, and extracellular matrix studies.

[Learn More](#)

Parameter	Value
Product Name	CTGF Rabbit pAb
Remarks/Alias	CTGF; NOV2; HCS24; IGFBP8
Species	Human
Gene ID	1490 (Human)
Immunogen	Synthetic peptide (aa 100-200 of human CTGF, NP_001892.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P29279
Protein Weight	38 kDa
Shipping	Ice bag