

Anti-Ctgf Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P29279

Item Number: CM0016908

Introduction

High-quality rabbit polyclonal antibody against human CTGF, validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, rat CTGF. Ideal for connective tissue research, chondrocyte differentiation, and osteogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	CTGF Rabbit pAb
Remarks/Alias	CTGF; NOV2; HCS24; IGFBP8
Species	Human
GeneID (Human)	1490
GeneID	1490
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CTGF (NP_001892.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29279
Protein Weight	38kDa
Shipping	Ice bag