

Anti-Ctip1/Bcl11A Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H165

Item Number: CM0012184



Introduction

Rabbit monoclonal antibody targeting human CTIP1/BCL11A (Q9H165). Validated for WB, IHC-P, ELISA and cross-reacts with rat. Ideal for research on globin switching, B-cell development, and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CTIP1/BCL11A Rabbit mAb
Remarks/Alias	EV19; CTIP1; DILOS; ZNF856; HBFQTL5; SMARCM1; 1A
Species	Human
GeneID (human)	53335
GeneID	53335
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 401-589 of human CTIP1/BCL11A (Q9H165)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H165
Protein Weight	91kDa
Shipping	Ice bag