

Anti-Ctip2/Bcl11B Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9C0K0

Item Number: CM0003336

Introduction

Recombinant rabbit monoclonal antibody against CTIP2/BCL11B, validated for Western blot, IHC-P, and ELISA. Reacts with human and mouse. Ideal for studying T-cell development and leukemia.

[Learn More](#)

Parameter	Value
Product Name	CTIP2/BCL11B Rabbit mAb
Remarks/Alias	ATL1; RIT1; CTIP2; IMD49; CTIP-2; IDDFSTA; SMARCM2; ZNF856B; ATL1-beta; ATL1-alpha; ATL1-delta; ATL1-gamma; hRIT1-alpha; CTIP2/BCL11B
Species	Human
GeneID (Human)	64919
GeneID	64919
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9C0K0
Protein Weight	96kDa
Shipping	Ice bag