

Anti-Ctip2/Bcl11B Rabbit Polyclonal Antibody For Wb - Q9C0K0

Item Number: CM0019304

Introduction

Rabbit polyclonal antibody targeting human CTIP2/BCL11B (Q9C0K0). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Key regulator of T-cell differentiation and survival, tumor suppressor, transcription repressor.

[Learn More](#)

Parameter	Value
Product Name	CTIP2/BCL11B Rabbit pAb
Remarks/Alias	ATL1; RIT1; CTIP2; IMD49; CTIP-2; IDDFSTA; SMARCM2; ZNF856B; ATL1-beta; ATL1-alpha; ATL1-delta; ATL1-gamma; hRIT1-alpha; CTIP2/BCL11B
Species	Human
Gene ID (Human)	64919
Gene ID	64919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-310 of human CTIP2/BCL11B (NP_612808.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0K0
Protein Weight	96kDa
Shipping	Ice bag