

Anti-Ctnna3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ui47

Item Number: CM0002155

Introduction

CTNNA3 rabbit monoclonal antibody validated for Western blot, immunohistochemistry, and ELISA in human, mouse, and rat samples. Suitable for studying cell adhesion, cardiomyopathy, and desmosome biology.

[Learn More](#)

Parameter	Value
Product Name	CTNNA3 Rabbit mAb
Remarks/Alias	VR22; ARVD13; CTNNA3
Species	Human
GeneID (Human)	29119
GeneID	29119
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-895 of human CTNNA3 (Q9UI47).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UI47
Protein Weight	100kDa
Shipping	Ice bag