

Anti-Ctnna3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Ui47

Item Number: CM0024645

Introduction

Rabbit polyclonal antibody against human CTNNA3 (Alpha T-catenin), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	CTNNA3 Rabbit pAb
Remarks/Alias	VR22; ARVD13; CTNNA3
Species	Human
GeneID (Human)	29119
GeneID	29119
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human CTNNA3 (NP_001278062.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UI47
Protein Weight	100kDa
Shipping	Ice bag