


Anti-Ctnnbip1 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Nsa3

Item Number: CM0003071



Introduction

Rabbit monoclonal antibody against human CTNNBIP1 (ICAT), a Wnt pathway negative regulator. Validated for WB and ELISA. Immunogen: N-terminal peptide (1-81); 9 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	CTNNBIP1 Rabbit mAb
Remarks/Alias	ICAT; CTNNBIP1
Species	Human
GeneID (Human)	56998
GeneID	56998
Immunogen	A synthetic peptide corresponding to amino acids 1-81 of human CTNNBIP1 (Q9NSA3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NSA3
Protein Weight	9kDa
Shipping	Ice bag