

Anti-Ctnnbl1 Polyclonal Antibody For Wb, If/Icc, And Elisa - Q8Wya6

Item Number: CM0030416



Introduction

CTNNBL1 Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat. Suitable for splicing and somatic hypermutation research.

[Learn More](#)

Parameter	Value
Product Name	CTNNBL1 Rabbit pAb
Remarks/Alias	NAP; P14L; IMD99; PP8304; C20orf33; dj633020.1; CTNNBL1
Species	Human
Gene ID	56259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 264-563 of human CTNNBL1 (NP_110517.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WYA6
Protein Weight	65kDa
Shipping	Ice bag